



Working Together to lighten the load

DOWNDRAFT AUTOPSY TABLE, MODEL RCS



- **STYLISH ATTRACTIVE FUNCTIONAL DESIGN**
- **SECTIONAL REVOMABLE WORKTOP**
- **POSITIVE AIRFLOW AROUND & OVER TABLE**
- **180° ROTATION WITH ELEC. HEIGHT ADJUSTMENT**
- **ELECTRICAL INTERLOCKED BRAKE MECHANISM**
- **EASY TO CLEAN, CONCEALED ACCESS TO PLENUM**
- **UNIQUE PATENTED 'ANTI WOBBLE' TABLE LOCK**
- **304 STAINLESS FABRICATION**

Featuring an aesthetically pleasing, slim pedestal, this stainless steel autopsy table is designed to provide perimeter downdraft ventilation, ensuring full performance when in use.

The table offers simple user electrical height adjustment, elevating the entire chassis a nominal 250mm, whilst the central pedestal accommodates the engineering services.

The worktop is interlocked and located onto the chassis and can be fabricated in either two or three sections to simplify removal and cleaning routines. The unique design also provides a positive table lock essential for safe operation, eliminating the 'wobble' effect found on traditional designs.

The fabrication is mounted on a baseplate, designed for fixing to the sub floor, through which the extract and services are delivered.

Easily adapted with many flexible features this unique, patented and bespoke design brings another innovative design to the mortuary, successfully integrating downdraft technology.



DOWNDRAFT AUTOPSY TABLE

OPTIONS

- BESPOKE WORKTOPS TO SUIT REQUIREMENTS
- EXTERNAL LP WATER REGULATOR
- THERMOSTATIC MIXER
- H&C SPRAY FACILITY
- PERIMETER BARS FOR INSTRUMENT CLAMPING
- 316 GRADE FABRICATION
- FULL DESIGN AND INSTALL

FEATURES

- Stylish attractive yet functional design
- Sectional removable worktop
- Perimeter extract with solid worktop
- Positive airflow around and over table
- 180° Rotation with electrical height adjustment
- Electrical interlocked brake mechanism
- Unique patented positive 'anti wobble' table lock
- Easy to clean, concealed access to plenum
- 304 grade stainless fabrication
- Internal sump under worktop fitted with sparge rinse facility

LEEC Ltd
Private Road No 7
Colwick Ind. Est.
Nottingham
NG4 2AJ

Tel: 0044 (0) 115 961 6222
Fax: 0044 (0) 115 961 6680

E-mail: sales@leec.co.uk
Website: www.leec.co.uk

REF
C19-03