

DISSECTING BENCH, MODEL VDB



- FULLY VENTILATED WORKTOP
- WIDE, ACCESSIBLE PLATFORM
- EASY TO CLEAN, NO HEAVY PARTS
- LEFT OR RIGHT HANDED LAYOUT
- STAINLESS STEEL WORKTOP
- 750MM DEEP, VARIOUS LENGTHS
- PLAIN OR PERFORATED WORKTOP

A stainless steel 'downdraft' ventilated bench, featuring a removable upper sectional worktop. Extract is provided through a central spigot within a sump situated under the worktop, a sink bowl can be sited at either end.

The unit is available with a number of options, suited typically for Laboratory or Mortuary use, but with other applications.

Based on a standard format but with a number of sizes available, these benches can be tailored to the clients exact needs.



FEATURES

- Plain or Perforated removable stainless work surface
- Extract upstand, with adjustable splash shield
- Downdraft extract with stainless spigot outlet for onward duct connection
- Fully profiled perimeter
- Various size sink bowls, right or left handed
- 150mm long. Lever action, biflo mixer taps
- Water sparge under worktop for continuous wash
- Fitted with adjustable feet and rear upstand
- Peninsular, wall or cloaked versions

OPTIONS

- Var. pressure trigger action spray
- Adjustable height
- Hand operated spray facility
- Integrated or separate sluice unit
- Inner removable organ bowl
- Pathological macerator
- Full ventilation design & installation
- Quick release water couplings
- Formalin Dispensing equipment
- Polypropylene cutting boards
- St st waste pipe work
- Airflow alarms
- Rear upstands
- Cloaking panels

DOWNDRAFT DISSECTING BENCH

REF C5-08

LEEC Ltd
 Private Road No 7
 Colwick Ind. Est.
 Nottingham
 NG4 2AJ

Tel: 0044 (0) 115 961 6222
 Fax: 0044 (0) 115 961 6680

E-mail: sales@leec.co.uk
 Website: www.leec.co.uk

